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# EN #291322 AniScale Animal Scale Service Manual

You can change basic settings and calibrate the animal scale in the setup menu. Only change settings after consultation with the manufacturer.

## 1. Operation

### 1.1 Button functions

Button	Function
SETUP	Open the setup menu (hold down when switching on)
ESC	Exit setup.
ENTER	Discard changes. Next menu item. Confirm permanent change
Arrow keys (▲ / ▼)	Change current parameters.

### 1.2 Open setup menu

The device is off. Hold down the setup button when switching the device on. The display will show *PASS*. You will now be asked for your password. Change the displayed value to '2' using the arrow keys and confirm with ENTER to open the setup menu. The display will show *SETUP*.

### 1.3 Changing parameters

- You can scroll through the parameters in the setup menu using ENTER (see table in paragraph 3). You can change each value using the arrow keys.
- If you move to the next parameter using ENTER, the changed value for this parameter will be permanently saved. The display will show *End*.
- If you do not want to make the change permanent, then discard the change with ESC. The display will show *Exit*.

### 1.4 Close setup menu

You can exit the setup menu using ESC

## 2. Calibrating the animal scale

Proceed as follow to calibrate the animal scale:

- Open the setup menu as described above.
- Set the amplification of the AD converter (Gain). The amplification should be selected so that the measured analogue value under full scale load is around 30,000 (the measured analogue value is also displayed in the lower display line in parallel to the set amplification). The value may not exceed 32,767. Confirm the entry with ENTER.
- Menu item CAL2 will be shown on the display. Press the ▲ button to start the calibration. The measured analogue value will be displayed. Once the value has stabilised, press the ENTER key
- You will be taken to menu point SCL2. Enter the weight you used for calibration under CAL2 here. Confirm the entry with the ENTER key.
- Menu item CAL1 will now be shown on the display. Remove the load from the scale. Press the ▲ button to start the calibration. The analogue value measured on the empty scale will be displayed. Once the value has stabilised, press the ENTER key to confirm.
- Press the ▲ button once in menu item SCL1; leave the value set to '0'. Confirm with the ENTER key.
- Set the resolution of the display in menu item *RESOL*. The resolution must be greater than 1/4000 of the measuring range.
- All other settings should also be sensibly specified. They do not, however, affect calibration.

If the calibration does not work, or the measured analogue value does not change with the weight, please contact the manufacturer.

## 3. Parameters in the SETUP menu

Anzeige	Einstellbereich	Parameter	Bedeutung
<i>GA in</i>	1 - 7	Gain	Set the amplification of the AD amplifier
<i>CAL2</i>	cannot be set	Calibration point 2	The calibration weight is on the scale
<i>SCL2</i>	0 - 999999	Calibration weight 2	Weight for the calibration on calibration point 2
<i>CAL1</i>	cannot be set	Calibration point 1	The calibration weight has been removed
<i>SCL1</i>	0	Calibration weight 1	Load weight at calibration point 1 (0 kg)
<i>RESOL</i>	1 - 250	Resol	Weight display resolution (default 1.2.5)
<i>dCPT</i>	0 - 0.000	Decimal pint	Setting the decimal point
<i>FLtr1</i>	0 - 99	Filter 1	Influences weight stabilisation when changing the weight (default 60 - 90)
<i>FLtr2</i>	0 - 250	Filter 2	Filters the changes around a stable weight (default approx. 20)
<i>AVrg</i>	0 - 50	Avrg	Number of measurements used to calculate an average
<i>our 9d</i>	0 - 999999	Overload	Weight limit for displaying the overload warning, but continued display of the measurement value
<i>our MAX</i>	0 - 999999	Overload Max	Weight limit for overload, display – in the display (default 110% of the capacity)
<i>MANCL2</i>	0 - 999999	Manual calibration point 2	Manual setting of the calibration point 2. Is used to override the calibration
<i>MANCL1</i>	0 - 999999	Manual calibration point 1	Manual setting of the calibration point 1. Is used to override the calibration
<i>MANSEt</i>	Yes - No	Factory settings	Reset to factory settings
<i>PEAK</i>	Yes - No	Peak	The second line shows the highest measured weight
<i>RS485</i>	Yes - No	RS485	Communication active via RS485 interface
<i>KEYb</i>	KeyB6 / KeyB8	KeyB	Use keyboard with 6 or 8 keys
<i>Print</i>	Yes - No	Print	Printer active
<i>Yoh</i>	00 - 99	Year	Current year
<i>MEs</i>	01 - 12	Month	Current month
<i>dEn</i>	1 - 31	Day	Current day
<i>hod</i>	0 - 23	Hour	Current time - hour
<i>Min</i>	0 - 59	Minute	Current time - minute
<i>SEc</i>	0 - 59	Second	Current time - second
<i>COUNT</i>	0 - 255	Number of measurements	Number of measurements used to calculate a display. Impacts the weighing time
<i>Min</i>	0 - 9999	Minimum weight	Anything under this weight will reset the scale and allow a new weighing process

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